

ELECTROMAGNETIC VISCOMETER

EV 1300

Precision viscosity measurement
under extreme pressure & temperature conditions



0.2 – 10,000 cP

Viscosity Range



Up to 200°C

Temperature



Up to 20,000 psi

Pressure



20 cc

Sample Volume

EV 1300 – Vibrating Rod Viscometer

Technology & Specifications

OPERATION

TECHNICAL SPECIFICATIONS

Viscosity Range	0.2 – 10,000 cP
Accuracy	±2% of displayed value
Pressure Range	Up to 20,000 psi
Temperature Range	Ambient to 200°C
Sample Volume	20 cc
Wetted Parts	Stainless steel
Power Supply	220 VAC
Dimensions (W×H×D)	700 × 950 × 650 mm



OPERATING PRINCIPLE

- ❑ A rod vibrates inside a 20-ml measurement cell at a known frequency.
 - ❑ A piezoelectric sensor measures this damping. The electrical signal is proportional to the fluid's viscosity
 - ❑ The vibration amplitude changes according to the viscosity of the surrounding fluid
 - ❑ The viscous fluid dampens the vibration. The higher the viscosity, the greater the damping.
 - ❑ The viscometer converts this signal into a viscosity value. Measurement is continuous and real-time.
-
- ❑ Two interchangeable probes cover the full range: Probe 1 for 0.2–200 cP, Probe 2 for 200–10,000 cP.
 - ❑ Inlet & outlet valves are provided within the cell
 - ❑ An integrated air bath maintains uniform thermal control up to 200°C.
 - ❑ Factory-calibrated; field recalibration possible using silicone oil (NIST standards)

