



Low pressure syringe pump



PV series



Features

- ❑ Single cylinder configuration
- ❑ Ambient temperature
- ❑ Stainless steel wetted parts
- ❑ Supervision software
- ❑ Ethernet connection for data logging
- ❑ Horizontal / vertical position

Large range of applications:

- ❑ Water pore volume change
- ❑ Pressure control
- ❑ Soil studies

Model	Pressure	Volume	Max flow rate	Mini Flow rate
	(psi)	(ml)	(ml/min)	(ml/min)
PV 200	600	200	70	0.0001
PV 500	225	500	150	0.0001